

HMG20A Antibody



PACO09742

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HMG20A

Uniprot

Q9NP66

Synonyms:

high-mobility group 20A; HMG20A; FLJ10739; HMGX1; HMGXB1 ;

Immunogen:

Human HMG20A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A