

HMGB4 Antibody



PACO09746

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HMGB4

Uniprot

Q8WW32

Synonyms:

high-mobility group box 4; HMGB4; FLJ40388; MGC88128; dJ1007G16.5 ;

Immunogen:

Human HMGB4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A