

HMGN1 Antibody



PACO09749

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HMGN1

Uniprot

P05114

Synonyms:

high-mobility group nucleosome binding domain 1; HMGN1; FLJ27265; FLJ31471; HMG14; MGC104230; MGC117425 ;

Immunogen:

Human HMGN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A