

# HMOX2 Antibody



PACO09754

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

HMOX2

**Uniprot**

P30519

**Synonyms:**

heme oxygenase (decycling) 2; HMOX2; HO-2 ;

**Immunogen:**

Human HMOX2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A