

PACO09761

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

HNRNPC

**Uniprot**

P07910

**Synonyms:**

heterogeneous nuclear ribonucleoprotein C (C1/C2); HNRNPC; C1; C2; HNRNP; HNRPC; MGC104306; MGC105117; MGC117353; MGC131677; SNRPC ;

**Immunogen:**

Human HNRNPC.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A