

HNRNPD Antibody



PACO09762

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HNRNPD

Uniprot

Q14103

Synonyms:

heterogeneous nuclear ribonucleoprotein D (AU-rich element RNA binding protein 1, 37kDa); HNRNPD; AUF1; AUF1A; HNRPD; P37; hnRNP0 ;

Immunogen:

Human HNRNPD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A