

HNRNPH1 Antibody



PACO09764

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HNRNPH1

Uniprot

P31943

Synonyms:

heterogeneous nuclear ribonucleoprotein H1 (H); HNRNPH1; DKFZp686A15170; HNRPH; HNRPH1; hnRNPH ;

Immunogen:

Human HNRNPH1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A