

# HNRNPUL1 Antibody



PACO09768

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HNRNPUL1

**Uniprot**

Q9BUJ2

**Synonyms:**

heterogeneous nuclear ribonucleoprotein U-like 1; HNRNPUL1; E1B-AP5; E1BAP5; FLJ12944; HNRPUL1 ;

**Immunogen:**

Human HNRNPUL1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A