

# HOXB1 Antibody



PACO09776

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HOXB1

**Uniprot**

P14653

**Synonyms:**

homeobox B1; HOXB1; HOX2; HOX2I; Hox-2.9; MGC116843; MGC116844; MGC116845  
;

**Immunogen:**

Human HOXB1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A