

HOXC10 Antibody



PACO09777

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HOXC10

Uniprot

Q9NYD6

Synonyms:

homeobox C10; HOXC10; HOX3l; MGC5259 ;

Immunogen:

Human HOXC10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A