

HOXD13 Antibody



PACO09780

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HOXD13

Uniprot

P35453

Synonyms:

homeobox D13; HOXD13; BDE; BDSD; HOX4l; SPD ;

Immunogen:

Human HOXD13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A