

HOXD3 Antibody



PACO09781

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

HOXD3

Uniprot

P31249

Synonyms:

homeobox D3; HOXD3; HOX1D; HOX4; HOX4A; Hox-4.1; MGC10470 ;

Immunogen:

Human HOXD3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A