

PACO09785

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HPDL

Uniprot

Q96IR7

Synonyms:

4-hydroxyphenylpyruvate dioxygenase-like; HPDL;4-HPPD-L; GLOXD1; MGC15668; RP4-534D1.1 ;

Immunogen:

Human HPDL.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A