

# HS3ST1 Antibody



PACO09802

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HS3ST1

**Uniprot**

O14792

**Synonyms:**

heparan sulfate (glucosamine) 3-O-sulfotransferase 1; HS3ST1;3OST;3OST1 ;

**Immunogen:**

Human HS3ST1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A