

# HSD11B1 Antibody



PACO09804

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HSD11B1

**Uniprot**

P28845

**Synonyms:**

hydroxysteroid (11-beta) dehydrogenase 1; HSD11B1;11-DH;11-beta-HSD1; HDL; HSD11; HSD11B; HSD11L; MGC13539; SDR26C1 ;

**Immunogen:**

Human HSD11B1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A