

HSD17B11 Antibody



PACO09807

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HSD17B11

Uniprot

Q8NBQ5

Synonyms:

hydroxysteroid (17-beta) dehydrogenase 11; HSD17B11;17-BETA-HSD11;17-BETA-HSDXI; DHRS8; PAN1B; RETSDR2; SDR16C2 ;

Immunogen:

Human HSD17B11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A