

HSD17B14 Antibody



PACO09808

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HSD17B14

Uniprot

Q9BPX1

Synonyms:

270aa,28kd;

Immunogen:

Human HSD17B14.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A