

HSF2 Antibody



PACO09818

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HSF2

Uniprot

Q03933

Synonyms:

heat shock transcription factor 2; HSF2; MGC117376; MGC156196; MGC75048 ;

Immunogen:

Human HSF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A