

HSF4 Antibody



PACO09819

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HSF4

Uniprot

Q9ULV5

Synonyms:

heat shock transcription factor 4; HSF4; CTM ;

Immunogen:

Human HSF4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A