

HSPA9 Antibody



PACO09829

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HSPA9

Uniprot

P38646

Synonyms:

heat shock 70kDa protein 9 (mortalin); HSPA9; CSA; GRP75; HSPA9B; MGC4500; MOT; MOT2; MTHSP75; PBP74; mot-2 ;

Immunogen:

Human HSPA9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A