

HSPD1 Antibody



PACO09835

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HSPD1

Uniprot

P10809

Synonyms:

heat shock 60kDa protein 1 (chaperonin); HSPD1; CPN60; GROEL; HLD4; HSP60; HSP65; HuCHA60; SPG13 ;

Immunogen:

Human HSPD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A