

HSPH1 Antibody



PACO09836

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

HSPH1

Uniprot

Q92598

Synonyms:

heat shock 105kDa/110kDa protein 1; HSPH1; DKFZp686M05240; HSP105; HSP105A; HSP105B; KIAA0201; NY-CO-25 ;

Immunogen:

Human HSPH1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A