

ICA1L Antibody



PACO09850

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ICA1L

Uniprot

Q8NDH6

Synonyms:

islet cell autoantigen 1,69kDa-like; ICA1L; ALS2CR14; ALS2CR15; DKFZp434E1919; MGC138440 ;

Immunogen:

Human ICA1L.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A