

IDI2 Antibody



PACO09867

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

IDI2

Uniprot

Q9BXS1

Synonyms:

isopentenyl-diphosphate delta isomerase 2; IDI2; IPP12 ;

Immunogen:

Human IDI2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A