

IDUA Antibody



PACO09869

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

IDUA

Uniprot

P35475

Synonyms:

iduronidase, alpha-L-; IDUA; IDA; MPS1 ;

Immunogen:

Human IDUA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A