

IFI30 Antibody



PACO09871

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IFI30

Uniprot

P13284

Synonyms:

interferon, gamma-inducible protein 30; IFI30; GILT; IFI-30; IP30; MGC32056 ;

Immunogen:

Human IFI30.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A