

IFIT3 Antibody



PACO09873

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

IFIT3

Uniprot

O14879

Synonyms:

interferon-induced protein with tetratricopeptide repeats 3; IFIT3; CIG-49; GARG-49; IFI60; IFIT4; IRG2; ISG60; RIG-G ;

Immunogen:

Human IFIT3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A