

IFRD2 Antibody



PACO09884

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IFRD2

Uniprot

Q12894

Synonyms:

interferon-related developmental regulator 2; IFRD2; FLJ40446; IFNRP; SKMc15; SM15 ;

Immunogen:

Human IFRD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A