

IFT52 Antibody



PACO09887

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

IFT52

Uniprot

Q9Y366

Synonyms:

intraflagellar transport 52 homolog (Chlamydomonas); IFT52; C20orf9; CGI-53; NGD5; dJ1028D15.1 ;

Immunogen:

Human IFT52.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A