

IGF2BP3 Antibody



PACO09893

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

IGF2BP3

Uniprot

O00425

Synonyms:

insulin-like growth factor 2 mRNA binding protein 3; IGF2BP3; DKFZp686F1078; IMP-3; IMP3; KOC1; VICKZ3 ;

Immunogen:

Human IGF2BP3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A