

# IGFBP2 Antibody



PACO09895

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IGFBP2

**Uniprot**

P18065

**Synonyms:**

insulin-like growth factor binding protein 2, 36kDa; IGFBP2; IBP2; IGF-BP53 ;

**Immunogen:**

Human IGFBP2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A