

IKZF1 Antibody



PACO09905

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IKZF1

Uniprot

Q13422

Synonyms:

IKAROS family zinc finger 1 (Ikaros); IKZF1; Hs.54452; IK1; IKAROS; LYF1; PRO0758; ZNFN1A1; hlk-1 ;

Immunogen:

Human IKZF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A