

IL10RB Antibody



PACO09909

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IL10RB

Uniprot

Q08334

Synonyms:

interleukin 10 receptor, beta; IL10RB; CDW210B; CRF2-4; CRFB4; D21S58; D21S66; IL-10R2 ;

Immunogen:

Human IL10RB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A