

IL15RA Antibody



PACO09914

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IL15RA

Uniprot

Q13261

Synonyms:

interleukin 15 receptor, alpha; IL15RA; MGC104179 ;

Immunogen:

Human IL15RA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A