

IL17RB Antibody



PACO09918

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

IL17RB

Uniprot

Q9NRM6

Synonyms:

interleukin 17 receptor B; IL17RB; CRL4; EVI27; IL17BR; IL17RH1; MGC5245 ;

Immunogen:

Human IL17RB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A