

IL17RC Antibody



PACO09919

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IL17RC

Uniprot

Q8NAC3

Synonyms:

interleukin 17 receptor C; IL17RC; FLJ95963; FLJ96005; IL17-RL; IL17RL; MGC10763 ;

Immunogen:

Human IL17RC.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A