

IL22RA1 Antibody



PACO09930

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IL22RA1

Uniprot

Q8N6P7

Synonyms:

interleukin 22 receptor, alpha 1; IL22RA1; CRF2-9; IL22R; IL22R1 ;

Immunogen:

Human IL22RA1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A