

IL2RB Antibody



PACO09932

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

IL2RB

Uniprot

P14784

Synonyms:

interleukin 2 receptor, beta; IL2RB; CD122; P70-75 ;

Immunogen:

Human IL2RB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A