

IMMP2L Antibody



PACO09945

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

IMMP2L

Uniprot

Q96T52

Synonyms:

IMP2 inner mitochondrial membrane peptidase-like (S. cerevisiae); IMMP2L; IMP2; IMP2-LIKE ;

Immunogen:

Human IMMP2L

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A