

IMPA1 Antibody



PACO09949

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

IMPA1

Uniprot

P29218

Synonyms:

inositol(myo)-1(or 4)-monophosphatase 1; IMPA1; IMPA ;

Immunogen:

Human IMPA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A