

IMPDH2 Antibody



PACO09950

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

IMPDH2

Uniprot

P12268

Synonyms:

IMP (inosine monophosphate) dehydrogenase 2; IMPDH2; IMPD2; IMPDH-II ;

Immunogen:

Human IMPDH2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A