

ING3 Antibody



PACO09954

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ING3

Uniprot

Q9NXR8

Synonyms:

inhibitor of growth family, member 3; ING3; Eaf4; FLJ20089; ING2; p47ING3 ;

Immunogen:

Human ING3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A