

INHBB Antibody



PACO09958

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

INHBB

Uniprot

P09529

Synonyms:

inhibin, beta B; INHBB; MGC157939 ;

Immunogen:

Human INHBB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A