

# INTS3 Antibody



PACO09972

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

INTS3

**Uniprot**

Q68E01

**Synonyms:**

integrator complex subunit 3; INTS3; C1orf193; C1orf60; DKFZp686E1950; DKFZp686O20115; DKFZp781I1253; FLJ21919; FLJ44849; INT3 ;

**Immunogen:**

Human INTS3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A