

PACO09978

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

INVS

**Uniprot**

Q9Y283

**Synonyms:**

inversin; INVS; INV; KIAA0573; MGC133080; MGC133081; NPH2; NPHP2 ;

**Immunogen:**

Human INVS.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A