

IPO13 Antibody



PACO09982

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

IPO13

Uniprot

O94829

Synonyms:

importin 13; IPO13; IMP13; KAP13; KIAA0724; LGL2; RANBP13 ;

Immunogen:

Human IPO13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A