

IRF8 Antibody



PACO09992

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

IRF8

Uniprot

Q02556

Synonyms:

interferon regulatory factor 8; IRF8; H-ICSBP; ICSBP; ICSBP1; IRF-8 ;

Immunogen:

Human IRF8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A