

ISL1 Antibody



PACO09998

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ISL1

Uniprot

P61371

Synonyms:

ISL LIM homeobox 1; ISL1; ISLET1; Isl-1 ;

Immunogen:

Human ISL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A