

ITGB1BP1 Antibody



PACO10007

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ITGB1BP1

Uniprot

O14713

Synonyms:

integrin beta 1 binding protein 1; ITGB1BP1; DKFZp686K08158; ICAP-1A; ICAP-1B; ICAP-1alpha; ICAP1; ICAP1A; ICAP1B ;

Immunogen:

Human ITGB1BP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A