

ITGB3 Antibody



PACO10009

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ITGB3

Uniprot

P05106

Synonyms:

integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61); ITGB3; CD61; GP3A; GPIIIa ;

Immunogen:

Human ITGB3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A